

## 离心机选型指南\_2021CN

### 1. 常温、冷冻与加热离心机

常温	SPROUT	54	55	56	BASIC	DIAGNOSTIC
	SCIENCE	UNIVERSAL	260	351e	352	380
冷冻	150R	260R	260RH	352R	352RH	380R
加热	260RH	352RH				

### 2. 转速和离心力

低速	4,000rpm	BASIC	351e			
	6,000rpm	54	56	DIAGNOSTIC	SPROUT	
高速	15,000rpm	55	150R			
	18,000rpm	SCIENCE	UNIVERSAL	260	260R	
		260RH	352	352R	352RH	
		380	380R			

### 3. 转子应用

角转子						
0.2ml	150R	SCIENCE	UNIVERSAL	260/R/RH	352/R/RH	380/R
1.5/2ml	55	56	150R	SCIENCE	UNIVERSAL	260/R/RH
	352/R/RH	380/R				
1.8ml	150R	SCIENCE	UNIVERSAL	260/R/RH		
2ml 滤柱	150R	SCIENCE	UNIVERSAL	352/R/RH	380/R	
5ml	150R	SCIENCE	UNIVERSAL	260/R/RH	352/R/RH	
10/15ml	54	56	150R	BASIC	DIAGNOSTIC	UNIVERSAL
	260/R/RH	351e	352/R/RH	380/R		
25/30ml	DIAGNOSTIC	UNIVERSAL	260/R/RH	351e	352/R/RH	380/R
50ml	DIAGNOSTIC	UNIVERSAL	UNIVERSAL	260/R/RH	352/R/RH	380/R
60ml	BASIC	DIAGNOSTIC				
85ml	352/R/RH	380/R				



100ml	UNIVERSAL	260/R/RH	352/R/RH			
250ml	380/R					
BABCOCK	352/R/RH					
<b>水平转子</b>						
10ml	DIAGNOSTIC	UNIVERSAL				
10/15ml	BASIC	DIAGNOSTIC	UNIVERSAL			
40ml	352/R/RH					
60ml	BASIC					
100ml	DIAGNOSTIC	UNIVERSAL	260/R/RH			
200ml	351e	352/R/RH				
250ml	352/R/RH					
500ml	380/R					
750ml	380/R					
ASTM	380/R					
<b>酶标板转子</b>						
2M	DIAGNOSTIC	UNIVERSAL	260/R/RH			
6M/2D	352/R/RH					
8M/2D	351e	352/R/RH	380/R			
8M/4D	380/R					
10M/4D	380/R					
<b>毛细管转子</b>						
50mm	55					
75mm	150R	SCIENCE	UNIVERSAL	260/R/RH	352/R/RH	
<b>生物安全转子 CYTO</b>						
CYTO	DIAGNOSTIC	352/R/RH				

HEAD OFFICE  
 9 PANTYGRAIGWEN ROAD  
 PONTYPRIDD MID GLAMORGAN  
 UNITED KINGDOM